

IMPORTANT: It is always recommended to follow cleaning instructions outlined in the user manual for each ZEISS instrument.



What solvents can be used to clean ZEISS device surfaces?

ZEISS recommends a disinfectant such as isopropyl alcohol. Isopropyl alcohol contains 70% alcohol, which is recommended by the Centers for Disease Control and Prevention (CDC). Please refer to the device user manual for a list of other appropriate cleaning and disinfecting agents.



Can the OCT lenses be cleaned and with what frequency?

Yes. As specified in section 11.3.4 of the CIRRUS user manual (Cleaning Heavy Contamination), lenses should be cleaned using only water, isopropyl alcohol and acetone. Up to 99% isopropyl alcohol can be used for cleaning CIRRUS optical surfaces.

Cleaning of optical surfaces can be done as much as needed (such as between patients), but there will be an increased risk of damage to the optical surfaces if hard particles are present on the surface and too much pressure is applied. Be sure to follow user manual cleaning protocols.




Can the HFA bowl be cleaned and how?

The cleaning of Humphrey® perimeters is unique because of the delicate surfaces of the perimeter's bowl and optics. We have detailed cleaning instructions on how to treat your Humphrey perimeter and the testing room, for your staff and your patient's safety. Also included is

new guidance on how
cleaned. However,
specifically, should the
clinic determine that
such cleaning is
desired.

[HFA Cleaning Guide](#) [\(/content/dam/med/ref_international/19_hfacleaning_guide.pdf\)](#) 

[Quickstart Guide: Cleaning your HFA perimeter](#) [\(/content/dam/med/ref_international/19_hfacleaning_quickstart.pdf\)](#) 

COVID-19 guidance: Cleaning your HFA perimeter



Quick links to cleaning instructions for ZEISS ophthalmic devi

OCT systems


[CIRRUS 6000 \(/content/dam/med/ref_international/corona/pdfs/ifu-cleaning-excerpt-cirrus-6000.pdf\)](/content/dam/med/ref_international/corona/pdfs/ifu-cleaning-excerpt-cirrus-6000.pdf)

[CIRRUS 5000/500 \(/content/dam/med/ref_international/corona/pdfs/ifu-cleaning-excerpt-cirrus-5000-500.pdf\)](/content/dam/med/ref_international/corona/pdfs/ifu-cleaning-excerpt-cirrus-5000-500.pdf)

[CIRRUS Photo \(/content/dam/med/ref_international/corona/pdfs/ifu-cleaning-excerpt_cirrus_photo.pdf\)](/content/dam/med/ref_international/corona/pdfs/ifu-cleaning-excerpt_cirrus_photo.pdf)

[PLEX](#)

[\(/content/dam/med/ref_international/corona/pdfs/cleaning_instructions_plex_elite.pdf\)](/content/dam/med/ref_international/corona/pdfs/cleaning_instructions_plex_elite.pdf)



Not all products, services or offers are approved or offered in every market and approved labeling and instructions may vary from one country to another. For country-specific product information, see the appropriate country website. Product specifications are subject to change in design and scope of delivery as a result of ongoing technical development.

 (<https://www.youtube.com/user/CarlZeissMeditec>)



(<https://www.instagram.com/zeissmedicaltechnology>)
(<http://www.zeiss.de/en>)





> [\(https://www.facebook.com/zeissmeditecglobal/\)](https://www.facebook.com/zeissmeditecglobal/)



[\(https://www.linkedin.com/company/carl-zeiss-meditec/\)](https://www.linkedin.com/company/carl-zeiss-meditec/)



[/zeissmeditec\)](https://twitter.com/zeissmeditec)

Further MED Websites

[Carl Zeiss Meditec AG](#)

[. \(https://www.zeiss.com/campaign/footer_tracking-ag-website_link\)](https://www.zeiss.com/campaign/footer_tracking-ag-website_link)

[Medical Technology for Patients](#)

[. \(https://www.zeiss.com/meditec/patient/int/home.html\)](https://www.zeiss.com/meditec/patient/int/home.html)

Customer Services

[Peer Insights](#)

[. \(/meditec/int/resource-center/peer_insights.html\)](https://www.zeiss.com/campaign/footer_tracking-insights.html)

[Online Shop](#)

[. \(https://www.zeiss.com/campaign/footer_tracking-shop-website_link\)](https://www.zeiss.com/campaign/footer_tracking-shop-website_link)

[E-Learning Center](#)

[. \(http://www.zeiss.com/campaign/footer_tracking-academy-website_link\)](http://www.zeiss.com/campaign/footer_tracking-academy-website_link)

[ZEISS Cataract Services](#)

[. \(https://www.zeiss.com/campaign/footer_tracking-cataract_service-website_link\)](https://www.zeiss.com/campaign/footer_tracking-cataract_service-website_link)

[ZEISS Newsletter](#)

[Registration](#)

[. \(https://pages.zeiss.com/Newsletter_subscription_center.html?utm_campaign=newsletter-subscription&utm_medium=int-homepage&utm_source=footer\)](https://pages.zeiss.com/Newsletter_subscription_center.html?utm_campaign=newsletter-subscription&utm_medium=int-homepage&utm_source=footer)

Further ZEISS Solutions

[ZEISS Vision](#)

[. \(https://www.zeiss.com/campaign/footer_tracking-vision-website_link\)](https://www.zeiss.com/campaign/footer_tracking-vision-website_link)

[ZEISS Microscopy](#)

[. \(https://www.zeiss.com/campaign/footer_tracking-microscopy-website_link\)](https://www.zeiss.com/campaign/footer_tracking-microscopy-website_link)

ZEISS Group

[Careers](#)

[International](#)

[. \(https://www.zeiss.com/campaign/careers-website_link\)](https://www.zeiss.com/campaign/careers-website_link)

[Press & Media](#)

[. \(https://www.zeiss.com/campaign/newsroom-website_link\)](https://www.zeiss.com/campaign/newsroom-website_link)

[Contact \(/meditec/int/about-us/contact-us/general-inquiry.html\)](/meditec/int/about-us/contact-us/general-inquiry.html)

| [Publisher \(/meditec/int/legal-information/company-information.html\)](/meditec/int/legal-information/company-information.html)

| [Legal Notice \(https://www.zeiss.com/meditec/legal-notice\)](https://www.zeiss.com/meditec/legal-notice)

[. \(http://www.zeiss.de/en\)](http://www.zeiss.de/en)

| [Data Protection \(https://www.zeiss.com/data-protection\)](https://www.zeiss.com/data-protection)





[Cookie Preferences](#)

[Back to MED Support Now\(/meditec/int/med-support-now.html\)](/meditec/int/med-support-now.html)

EN 38 173 0005III 